



The Model 2470 Shoulder CPM™

Technical Specifications

Adduction: 20 degrees
Abduction: 180 degrees

Internal/External Rotation: 0 to 110 degrees

Flexion: 180 degrees
Extension: 20 degrees

Synchronized Abduction and Rotation
(20—180 degrees and 0—110 degrees)

Speed: 1 to 3 degrees per second



The Model 2470 Shoulder CPM offers patients the safe, anatomical motion they need. The flexible range of motion options and patient lock-out feature prevent inadvertent changes to settings and ensure patient safety.

Designed to deliver superior comfort and versatile configurations, the Model 2470 is ideal for clinic and home shoulder rehabilitation.

Product Features

- Prevents stiffness by reducing edema and muscle contracture.
- Set up is simple and easy to understand by both patient and clinician.
- Manual jog mode allows set-up according to patient tolerance.
- The Model 2470 hand-held control pendant features a bold LED digital display and push button technology—easily understood by hospital staff and patients.
- Universal left/right design. Unit can be set up easily and quickly for conversion without the need for additional accessories.
- Patient compliance meter records the patient's treatment time.
- Modular components make transportation and shipping convenient.
- Shoulder CPM softgoods are constructed with the highest quality materials and provide a superb fit for the patient. Featuring deep-pile sheepskin and high-tensile strength Velcro brand strapping, Furniss softgoods offer superior comfort and support.
- Two-year manufacturer warranty.